

Minitron

An all-round genius in a small space. From standard experiments with microorganisms to complex cultivations of animal and plant cells, the Minitron has been designed to perform a comprehensive range of applications. Each Minitron is configured to your exact needs whilst the design ensures for a straightforward upgrade if required.

Dimensions and Weights



| Exterior dimensions | | | |
|---------------------------------------|-----------------------------|--|--|
| Height single unit (with rubber feet) | 700 mm | | |
| Height single unit (with base) | 812 mm | | |
| Height two units stacked (with base) | 1495 mm | | |
| Height base | 130 mm | | |
| Height rubber feet | 17 mm | | |
| Depth (with door closed) | 625 mm | | |
| Depth (with door open) | 1389 mm | | |
| Various | | | |
| Interior dimension (w x d x h) | 570 mm x 528 mm x 508 mm | | |
| Tray size | N (480 mm x 420 mm) | | |
| Additional space for ventilation | Back: 100 mm Side: 80 mm | | |

| Weight single unit without base frame and options | | |
|---|--------|--|
| Single unit 25 mm throw75 kg | | |
| Single unit 50 mm throw | 77 kg | |
| Weight options and accessories | | |
| Base frame | 7.5 kg | |
| Cooling unit | 9 kg | |
| Hygienic Humidity Control (HHC) | 3 kg | |
| CO ₂ control | 0.5 kg | |
| Analog output | 1 kg | |
| Universal tray | 2.5 kg | |



Shaker Drive / Rotation Speed

| Direction of rotation | Clockwise |
|---|---|
| Throw | 25 mm or 50 mm |
| Setting range (25 mm throw) | 20 min ⁻¹ to 400 min ⁻¹ |
| Setting range (50 mm throw) | 20 min ⁻¹ to 350 min ⁻¹ |
| Increment setpoint | 1 min ⁻¹ |
| Accuracy control (at maximum rota- tion speed, full scale) | ±1% |

Max. Rotation Speeds

| Single unit | 25 mm throw | 50 mm throw |
|-------------------|-----------------------|-----------------------|
| | 400 min ⁻¹ | 350 min-1 |
| Two units stacked | 25 mm throw | 50 mm throw |
| Top unit | 400 min ⁻¹ | 300 min ⁻¹ |
| Bottom unit | 400 min ⁻¹ | 350 min ⁻¹ |

Hygienic Humidity Control (HHC) (Option)

| General | | |
|-----------------------------|--------------|--|
| Setting range | 20 % to 85 % | |
| Increment setpoint | 1% | |
| Accuracy control | ± 3 % | |
| Water consumption (typical) | 5 g/h | |
| Max. temperature for use | 40 °C | |

| Reachable values | AT ¹⁾ | IT 2) | |
|-----------------------------------|------------------|-------|--------|
| Max. value (without condensation) | 20 °C | 37 °C | 75 %rH |
| | 25 °C | 37 °C | 75 %rH |
| Min. value (dehumidification) | 25 °C | 30 °C | 70 %rH |
| | 25 °C | 40 °C | 50 %rH |

¹⁾ AT = ambient temperature

²⁾ *IT = temperature in incubation chamber*

Temperature Control

| Setting range | 4 °C to 65 °C |
|---|---------------|
| Increment setpoint | 0.1 °C |
| Accuracy control 4 °C to 50 °C | ± 0.3 °C |
| Accuracy control > 50 °C | ± 0.5 °C |
| Temperature distribution, deviation: max. – min. ¹⁾ | ± 0.8 °C |
| Temperature distribution, max. deviation to value on display ¹⁾ | ± 0.8 °C |

¹⁾ at 37 °C in incubation chamber, on tray with 5 flasks

Lowest Attainable Temperature

| Configuration | Lowest attainable temp. |
|-----------------------------|---------------------------|
| Single unit without cooling | 5 °C above ambient temp. |
| Single unit with cooling | 16 °C below ambient temp. |

CO₂ Control (Option)

| Setting range | | 0.1 % to 20 % |
|---|------------------------------|---------------|
| Increment setpoint | | 0.1 % |
| Accuracy control ¹⁾ | 0 % to 5 % CO ₂ | 0.5 % |
| | 5 % to 10 % CO ₂ | 0.6 % |
| | 10 % to 15 % CO ₂ | 0.7 % |
| | 15 % to 20 % CO ₂ | 0.8 % |
| Gas consumption at 5 % CO ₂ (air vent open) | | 2 NL/h |
| Max. temperature for use | | 60 °C |

¹⁾ at 1013 hPa, 20 °C to 40 °C

Operating Conditions

| Load | | |
|------------------------------------|-----------------------------|--|
| Load max. | 12 kg | |
| Load optimal | 9.5 kg | |
| Ambient conditions | | |
| Ambient temperature | 10 °C to 32 °C | |
| Ambient humidity | 10 % to 85 % | |
| Altitude operating location | max. 2000 m above sea level | |
| Pollution degree as per EN 61010-1 | 2 | |
| Minimum distance side and back | 80 mm | |



Materials

| Housing | Polyurethane (PUR-IHS) with flame retardant |
|---------------------------------|---|
| Door | PUR-IHS, safety glass |
| Cover plate temperature control | Stainless steel (AISI 304) |
| Shaking table | Aluminium, anodised |
| Universal tray | Aluminium, anodised |

Various

| IP rating | | IP20 | |
|----------------|---------------|--------------------------|--|
| Sound pressure | | 35 dB(C) | |
| compressor | Version 230 V | R1234yf (quantity: 55 g) | |
| | Version 115 V | R134a | |

Interfaces

| Alarm connection | Stereo jack, 3 5 mm, allows to send system alarms to an external system. | |
|--------------------------------|--|--|
| Ethernet interface | RJ45, 10/100 Mbps Ethernet | |
| Analog output (optional) | 8 channels, 4 mA to 20 mA; allows to control the device and record data. | |
| Profibus DP gateway (optional) | Allows to connect the device to a SCADA system to control the device and record data. | |
| Modbus TCP gateway (optional) | | |

Electrical Connection and Power Values

| General | 230 V | 115 V |
|---|----------------------------|-------|
| Mains voltage | 230 V (±10 %) 115 V (±10 % | |
| Mains frequency | 50/60 Hz 60 Hz | |
| Max. power consumption (base unit) | 500 W | |
| Max. power consumption (all options) | 650 W | |
| Max. current consumption (all options) | 2.8 A | 5.6 A |
| Fuse (two 5 x 20 mm fuses, time lag) | 6.3 A | |



Connections/Utilities



| Pos. | Connection | Size | Pressure | Requirements |
|------|------------------------|---|--------------------|---|
| 1 | Demineralised water In | UNF 1/4-28 for hoses 1/8" (= 3.2 mm) | max. 2.0 bar | Water hardness (CaCO₃ equivalent): < 0.01 mmol L⁻¹ Dissolved solids: < 10 mg L⁻¹ Recommendation: Reverse osmosis water with a conductivity of approx. 5 µS cm⁻¹ or ultra-pure water/WFl. Do not use tap water, not even as an additive to ultra-pure water. |
| 2 | Discharge outlet | Internal thread G1/4", for hose Ø = 10 mm | N/A | N/A |
| 3 | CO ₂ In | Hose nozzle DN04, for hose Ø = 3 mm to 4 mm | 0.4 bar to 0.6 bar | For best efficiency, it is recommended to use a gas with a high CO ₂ concentration (e.g. 99.5 %). |

eve®



eve[®] is a platform software for planning, execution and analysis of bioprocesses. eve[®] allows you to record bioprocess data and store it in a central database. The software offers workflows from simple bioprocesses to the planning and execution of complex strategies with various phases.

eve[®] makes it possible to generate and store bioprocess knowledge. Various libraries for storing information on organisms and culture media are available. Thanks to soft-sensors, additional knowledge can be generated.

In addition to INFORS HT products, biotech machines and analysis devices from third-part manufacturers can be connected. This makes it possible to holistically control, monitor and analyse bioprocesses using a single software.

eve^{*} is installed on a centralised server. Access takes place via a browser, no client side installation is required. Bioprocess data is therefore available directly via the browser and independent of the operating system.

Various packages of the software are available. This makes it possible to adapt it to the individual needs and requirements of its users. eve[®] (in the premium version) is also suitable for working in a validated environment as per FDA CFR 21 Part 11.